

Accessories

- 3mm thick non-slip embossed aluminum plate improves the compressive and wear resistance of the cold room floor and extends its service life.
- PVC curtain can reduce the leakage of cold air when opening the door and reduce temperature fluctuation.

Parameter--Cold Room

Model	Monobloc						Split					
	WIC (10 cbm)	WIF (20 cbm)	WIC (30 cbm)	WIC (40 cbm)	Combi. 40 cbm WIC/WIF		WIF (20 cbm)	WIC (30 cbm)	WIC (40 cbm)	Combi. 40 cbm WIC/WIF		
					WIC (10 cbm)	WIC (10 cbm)				WIC (25 cbm)	WIF (15 cbm)	
Capacity (m³)	10	20	30	40	25	15	20	30	40	25	15	
WHO PQS Code	E001/003											
Refrigerant	R448A											
Defrost Mod	Hot gas defrost						Electrical heating					
Internal Temperature Range (°C)	2-8	-20	2-8	2-8	2-8	-20	-20	2-8	2-8	2-8	-20	
Power Supply (V/ph-Hz)	220/1N-50	380/3N-50	380/3N-50	380/3N-50	380/3N-50	380/3N-50	380/3N-50	380/3N-50	380/3N-50	380/3N-50	380/3N-50	
Power (W)/Unit	894	1750	1590	1590	1590	1750	1750	1590	1860	1590	1750	
Refrigeration Output (W) Unit	1425	2580	3700	3700	3700	2580	2580	3700	4600	3700	2580	
Condensation Temperature (°C)	43											
Density (Kg/Cbm)	40+/-2											
U Value (W/m²K)	0.17											
Insulation Thickness (mm)	100	120	100	100	120	120	120	100	100	120	120	

*Haier Biomedical reserves the right to change products and specifications without prior notice.

Walk-in Cold/ Freezer Room

Application Scenario

The storage temperature of 2-8°C and -20°C is suitable for vaccine storage in CDC, hospitals, inoculation sites and other places.



Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International

Solution of Vaccine Cold Room—— Mono Block Unit Cold Room



The inner temperature is kept constantly at 2 ~ 8°C (or -15 ~ -25°C).



Dual refrigeration systems, one system can backup the other system, extremely safe.



The refrigeration system has a variety of protection functions for start delay, overheat, overload, high and low pressure, over-voltage, under-voltage, etc.




Automatic defrosting by electric or hot gas.



Over temperature alarm, and power failure alarm.



PQS certified by the WHO.



SPECIFICATIONS:



Insulation System

- Stainless steel eccentric locking device is embedded under the insulation plate. The room is connected in hexahedral locking method which makes the room solid with good insulation effect
- Double-sided color steel plate, the intermediate insulation material is rigid high pressure foaming fluorine-free polyurethane. Stainless steel panel is optional
- Cold room plate thickness: 100mm; Freezer room plate thickness: 120mm; Steel plate thickness: 0.5mm
Polyurethane density: 42kg/m³; Thermal conductivity: ≤0.021W/(m.k); Adhesive strength: 0.14MPa; Compressive strength: 258kPa; Bending bearing capacity: 5.67mm
- Door dimensions 800*1975*100mm. Semi-buried design provides good insulation effect and the size of the door can be customized
- The door is equipped with anti-lock device to prevent the personnel from being locked in the warehouse

Refrigeration System

- Dual refrigeration system is safe and reliable, which will automatically switch to the other system when failure occurs to one
- Dual system will switch regularly so that uniform wear of refrigeration unit will prolong the service life of the unit
- Fluorine-free refrigerant: R448A, ODP: 0
- Automatic defrosting by electric or hot gas
- The technology of liquid-subcooling and exhaust-precooling via defrosting water is adopted to save energy

Shelves

- The main body of the shelf is designed with aluminum profile, and the partition is made of stainless steel. Overall strong and durable to prevent corrosion.
- The shelf has 5 layers of partitions, and the distance between the partitions can be adjusted



Illumination



- Equipped with three-proof LED lights, water-proof grade ≥ IP68
- The illumination on the vertical surface of the lowest shelf ≥ 150Lux