

Biosafety Screening Swab Booth

— **HCYC-1200**



Positive Pressure Protection Virus Isolation to Create a Safe Sample Swabbing Environment

Virus Isolation Inner Compartment

Positive pressure chamber, adjustable pressure, UV elimination, isolation operation.

Medical Fresh Air System

G4 pre-filter combined with U15 high efficiency filter to achieve medical filtration.

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
Tel: +86-0532-88935955
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



@haierbiomedicalint



Haier Biomedical
International

Scope of Application

Applicable for the safe sample swabbing of biological samples in hospitals, testing institutions, high-speed trains, railway stations, airports and highway junctions.

Product Advantages

Medical Fresh Air System

- Positive pressure in the booth, the pressure is adjustable, with pressure failure alarm, isolate viruses safely.

Air Filtration System

- Adopting G4 pre-filter and U15 high-efficiency filter, the filter element is easy to replace.

Virus Killing System

- Regular disinfection with ultraviolet remote control, convenient and worry-free self-service disinfection of intelligent induction glove disinfection machine, simple and efficient

Language Education Communication System

- Intercom device ensures convenient communication between medical staff and sampling staff



Ergonomic Design



Translucent Design

Extra large PMMA window is installed on the cabin sheet metal frame, and the booth is sealed and translucent



Lighting Design

LED light equipped inside for easy lighting



Caster Design

Four-corner roller design at the bottom of the booth, easy to move



Environmental Design

- The roof awning design can insulate heat and protect from the rain.
- Equipped with air conditioning, the temperature is adjustable, comfortable and pleasant.

Specifications

Model	Voltage (V/Hz)	Power (W)	Exterior Dimension (W*D*H mm)	Interior Dimension (W*D*H mm)	HEPA Class	Booth Sealing Strip
HCYC-1200	220/50	1000	1200*1400*2700	1190*1390*2100	G4+U15	German RAMPF raw material