

Haier

HUMANIZED DESIGN ACCURATE MEASUREMENT

MODEL
U80H



ELECTRONIC BLOOD PRESSURE MONITOR

HUMANIZED DESIGN HUMANIZED STYLE


Adopt Nordic style design
combine home design style with simple and
practical functional characteristics



LARGE SCREEN DISPLAY EASY TO READ

The screen of sphygmomanometer reaches 3.5 inches
and
the large font design is convenient for users to read

3.5
INCHES



SCIENTIFIC ALGORITHM ACCURATE MEASUREMENT

The built-in MCU chip bring accurate results
and comfortable experience

EASY TO USE EASY TO MEASURE

Humanized design ensures
that the measurement
process is simple.

the user is easy to use,
and easy to read
the measurement results.



ARRAY MEMORY EASY TO QUERY

It supports array memory,
and stores two user measurement data
into different memory groups,
which is convenient for users to query

ITEM

Product Size:103*103*65mm

Product Weight:230g

Memory:2*90

Measurement range

0-299 mmHg(0-37.30kPa)

Pulse:40-199 beats/min

Accuracy:Pressure ± 3 mmHg(0.4kPa) Pulse $\pm 5\%$ of the reading

Automatic power-off:180s