

Intelligent Vaccination Shelter

**Secure, Fast & Reliable Vaccination Facility
Guaranteed Vaccination Safety**



- Mobile vaccination shelter
- Rapid deployment and ready-to-use
- Clean and comfortable environment
- Long holdover time after power off, safe vaccine storage

Innovative & Ergonomic Design

- Transportable, can be deployed rapidly
- Intelligent monitoring system
- Multiple medical storage cabinets
- Customized layout to suit users needs
- Air disinfection system
- Shelter footprint is small

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
Tel: +86-0532-88935593
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International

Scope of Application

It is suitable for vaccination in vaccination centers, hospitals, CDC, schools, large enterprises and other institutions; the shelter is movable and applicable for vaccinating in large batches.

Product Advantages



Shelter Structure

The unit is low cost and mobile, can be deployed rapidly and it has small footprint.



Vaccine Storage Equipment

Equipped with an 8-small-drawers vaccination refrigerator with superior temperature uniformity and a solar direct drive refrigerator with long holdover time.



Storage Cabinet

Equipped with a number of storage cabinets for storage of inoculation cotton swabs, disposable needles and other items.



Intelligent Monitoring System

Independent monitoring of each operation area to achieve real-time supervision and inspection of the whole process of vaccination.



Air Disinfection System

Equipped with air purification sterilizer and ultraviolet disinfection, the cabin environment is safer.



Customized Internal Layout Scheme

Designed to suit multiple shelter schemes and can be customized according to user needs.

Top Layout View

The shelter is divided into three areas: registration area, vaccination area and cold chain area. Pre-check and registration of vaccination are conducted in the registration area. Vaccination is carried out in the vaccination area, and solar direct drive refrigerator is placed in the cold chain area.



Vertical View



Right View



Left View

Parameters

| Product Name | Power Supply (V/Hz) | Power (W) | Exterior Dimensions (W*D*H mm) | Interior Dimensions (W*D*H mm) | Weight (Kg) |
|---------------------------------|---------------------|-----------|--------------------------------|--------------------------------|-------------|
| Intelligent Vaccination Shelter | 220 (110)/50 (60) | 2,500 | 6,100*2,400*2,600 | 5,700*2,300*2,200 | 3,500 |