



Blood Bank Refrigerators

HXC-158

Scope of Application :

Applicable for the storage of blood and blood related products and other items that need to be stored at 2-6°C in hospitals and blood stations.

Innovative Design

- Frequency Conversion Compressor
- Adjustable Shelf And Basket
- Condensation Free Design
- Intelligent Alarm

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
Tel: +86-0532-88935593
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com



Product Advantages



Adjustable Shelves

The shelves are adjustable to any height to fulfill different requirements



Five Alarm Functions

High/low temperature, power failure, sensor error, low battery, door ajar



Three Types of Alarm

Audible buzzer, visible flashing light, remote alarm



Superior Performance

Adopting air cooling technology for the inside temperature uniformity and stability

Specifications



Model	HXC-158			
Technical Data	Type	Basket-Type		
	Climate class	ST		
	Cooling type	Forced air cooling		
	Defrost mode	Auto		
	Refrigerant	HC		
	Noise (dB(A))	42		
Performance	Temperature range (°C)	4±1		
	Ambient temperature (°C)	10~38		
Control	Controller	Microprocessor		
	Display	LED		
Electrical Data	Power supply (V/Hz)	220-240/50/60		
	Power (W)	320		
	Electrical current (A)	2.6		
Dimensions	Capacity (L/Cu.Ft)	158/5.6		
	Blood storage capacity (450ml blood bags)	84		
	Net/Gross weight (approx)	kg	107/120	
		lbs	235.9/264.6	
	Interior dimensions (W*D*H)	mm	460*370*950	
		in	18.1*14.6*37.4	
	Exterior dimensions (W*D*H)	mm	560*570*1530	
		in	22.0*22.4*60.2	
	Packing dimensions (W*D*H)	mm	645*675*1680	
		in	25.4*26.6*66.1	
Container load (20'/40'/40'H)	27/54/54			
Alarms	High/Low temperature	Y		
	Power failure	Y		
	Sensor error	Y		
	Low battery	Y		
	Door ajar	Y		
	Remote alarm	Y		
Accessories	Caster	Y		
	Foot	Y		
	Test hole	Y		
	Shelves/Basket	4/4		
	Drawers/Inner doors	-/2		
	USB interface	Optional		
	Temp recorder	Optional		
Others	Certification	CE		