

# MOBILE CLINIC with CT SCANNER UNIT

For imaging diagnosis of  
various systems of the whole  
body



- **ASSURE:** Medical electrical appliances are inspected by Haier Biomedical
- **SPECIAL CUSTOMIZATION:** Customized Appearance Painting Service
- **SEPARATE SPACE:** avoid cross infection
- **ACCEPTABLE CUSTOMIZATION:** equipment can be adjusted, purpose can be adjusted, design can be adjusted



扫描二维码关注海尔生物医疗  
官方网站, 让生命更美好



扫描二维码关注海尔生物医疗  
微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.  
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.  
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China  
URL: [www.haierbiomedical.com](http://www.haierbiomedical.com)  
Global unified hotline: +86-532-88935593

**Haier Biomedical**  
Make Life Better

### Features

The portable CT unit can be used for imaging diagnosis of diseases of various systems in the whole body, especially for lung examination during the epidemic, and it is easy to move and transport, reducing the risk of infection to the hospital for examination.



### Parameters

Overall dimensions of the vehicle (mm)	11985*2550*3980
Vehicle weight (kg)	30400
Chassis model	Shandeka ZZ5326XXYN466GF1K
engine model	MC07.31-60
Emission Standards	Chinese Standard VI
suspension system	air suspension
Cab occupants	3
Maximum speed (km/h)	89



扫描二维码关注海尔生物医疗  
官方网站, 让生命更美好



扫描二维码关注海尔生物医疗  
微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.  
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.  
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China  
URL: [www.haierbiomedical.com](http://www.haierbiomedical.com)  
Global unified hotline: +86-532-88935593

Haier Biomedical  
Make Life Better