

# MOBILE CLINIC GENERAL SURGERY UNIT

Equipped with vehicle-mounted generator set, equipped with independent working ability



- **ASSURE:** Medical electrical appliances are inspected by Haier Biomedical
- **SPECIAL CUSTOMIZATION:** Customized Appearance Painting Service
- **SEPARATE SPACE:** avoid cross infection
- **ACCEPTABLE CUSTOMIZATION:** equipment can be adjusted, purpose can be adjusted, design can be adjusted



扫描二维码关注海尔生物医疗  
官方网站, 让生命更美好



扫描二维码关注海尔生物医疗  
微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.  
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.  
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China  
URL: [www.haierbiomedical.com](http://www.haierbiomedical.com)  
Global unified hotline: +86-532-88935593

**Haier Biomedical**  
Make Life Better

### Features

As the core part of the mobile hospital, the general surgery unit is the operating room for medical rescue, equipped with a comprehensive operating bed, shadowless lamp and related equipment. Auxiliary devices such as oxygen supply system and air laminar flow purification device were installed, and the operating room reached ten thousand purification standards. Can realize the appendix excision, stomach excision, intestinal incision, thyroidectomy, hernia incision, cholecystectomy, tumor removal, fracture internal fixation surgery and gynecological delivery and other aseptic surgery.



### Parameters

Overall dimensions of the vehicle (mm)	11990*2500*3970
Vehicle weight (kg)	13500
Inner size (mm)	9307*2365*2105
Emission Standards	Chinese Standard V、GB3847-2005
Expand the box	Electronically controlled automatic expansion
Extended box width (mm)	900
Maximum speed (km/h)	89



扫描二维码关注海尔生物医疗官方网站, 让生命更美好



扫描二维码关注海尔生物医疗微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.  
 Hailemiao (Qingdao) Intelligent IoT Co., Ltd.  
 No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China  
 URL: [www.haierbiomedical.com](http://www.haierbiomedical.com)  
 Global unified hotline: +86-532-88935593