

MOBILE CLINIC - Ophthalmic Operating

Provide convenience for
cataract patients in remote
mountainous areas.



- **ASSURE:** Medical electrical appliances are inspected by Haier Biomedical
- **SPECIAL CUSTOMIZATION:** Customized Appearance Painting Service
- **SEPARATE SPACE:** avoid cross infection
- **ACCEPTABLE CUSTOMIZATION:** equipment can be adjusted, purpose can be adjusted, design can be adjusted



扫描二维码关注海尔生物医疗
官方网站, 让生命更美好



扫描二维码关注海尔生物医疗
微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.
Hailemiao (Qingdao) Intelligent IoT Co., Ltd.
No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China
URL: www.haierbiomedical.com
Global unified hotline: +86-532-88935593

Haier Biomedical
Make Life Better

Features

Haier Biomedical Eye Surgery Cart is a large-scale medical cart that integrates an operating room, a disinfection and cleaning room and an examination room. The car is equipped with professional inspection and screening equipment for eye disease screening. Let cataract patients enjoy professional medical services at their doorstep, and provide convenience for cataract patients in remote mountainous areas.



Parameters

Overall dimensions of the vehicle (mm)	9940*2470*3780
Vehicle weight (kg)	13500
Drive form	4X2
Chassis engine maximum power (kW)	155
Extended box width (mm)	650
generator set (kW)	6
Maximum speed (km/h)	89



扫描二维码关注海尔生物医疗官方网站, 让生命更美好



扫描二维码关注海尔生物医疗微信公众号, 让生命更美好

Qingdao Haier Biomedical Co., Ltd.
 Hailemiao (Qingdao) Intelligent IoT Co., Ltd.
 No. 280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China
 URL: www.haierbiomedical.com
 Global unified hotline: +86-532-88935593

Haier Biomedical
 Make Life Better