






Scope of Application

Used for smart storage for pick up of medications in hospital pharmacies. It can be connected with the HIS (hospital information system) so that pharmacists can manage dispensing corresponding medications according to physicians order.

Product Advantages

- 
High Efficiency
 2 separated conveyor systems, which are high efficiency
- 
Smart Indicating
 The location of medication will be reminded by LED indicators
- 
Grating Protection
 The equipment will stop once it detects any foreign bodies by optical grating
- 
Confirmation after Taking
 After taking the current drug, press "Confirm", the device will handle the next required drug
- 
Fast Dispense
 This device can handle multiple prescriptions so as to dispense fast and with efficiency



*The picture is for reference only, and the specific plan shall prevail.

Specifications

Model	Storage Location	Quantity of Drug Storage	Storage Method	Protection Measures	Power Supply (V/Hz)
HOH-ZN-108	108	More than 6000 boxes	The store closed, there is a sliding door at the window	Infrared light curtain protection	220-240~/50/60

*Haier Biomedical reserves the right to change products and specifications without prior notice.