



HYC-68/68A

Pharmacy Refrigerators



Audible and Visual Alarms

Audible buzzer and visible flashing light



Multiple Alarm Functions

High/low temperature, power failure, door ajar, sensor error, low battery, and remote alarm



Air Cooling, Auto Defrost

Superior cabinet temperature uniformity, auto defrost



Accurate Temperature Control



Interior LED



Lockable Design

Scope of Application

The under-counter pharmacy refrigerator offers ultra-low noise system, safety lock, complete alarm system for sample safety and protection. Widely used for storage of biological products, vaccines, drugs, reagents in pharmacies, pharmaceutical factories, hospitals, CDC, laboratories as key examples.

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International

Specifications



Technical Data

Cabinet Type	Upright
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	HC
Noise (dB (A))	41

Performance

Temp Range (C)	2~8
----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply (V/Hz)	220~240/50	115/60
Electrical Current (A)	0.85	1.7

Dimensions

Capacity (L/CuFt)	68/2.4
Net/Gross Weight (approx)	38/40 (kg) 83.8/88.2 (lbs)
Interior Dimensions (W*D*H)	415*385*505 (mm) 16.3*15.2*19.9 (in)
Exterior Dimensions (W*D*H)	495*580*660 (mm) 19.5*22.8*26.0 (in)
Packing Dimensions (W*D*H)	570*670*750(mm) 22.5*26.4*29.5(mm)
Container Load (20'/40'/40'H)	102/207/207



HYC-68



HYC-68A

Alarms

High/Low Temp	Y
Remote Alarm	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	N
Foot	Y/2
Porthole	Y
Shelves /Drawers	2/1
USB interface	Optional

Other

Certification	CE	UL
---------------	----	----

*Haier Biomedical reserves the right to change products and specifications without prior notice.