



**HYC-940F**

# Pharmacy Refrigerator



## Six Alarm Functions

High/low temperature, remote alarm, power failure, sensor error, low battery, door ajar



## Air Cooling, Auto Defrost

Better cabinet temperature uniformity, auto defrost



## Adjustable Shelves

The shelves are adjustable to any height to fulfill different requirements



## Self Closing Door

Energy conserving, maintains temperature uniformity to protect valuable samples

## Ergonomic Design



High Accuracy  
Temperature Control



Interior  
LED



Lockable  
Design

**Qingdao Haier Biomedical Co.,Ltd.**

No.280 Feng Yuan Road, High-tech Zone,  
Qingdao, 266109, P.R. China  
E-mail: [inquiry@haierbiomedical.com](mailto:inquiry@haierbiomedical.com)  
Website: [www.haiermedical.com](http://www.haiermedical.com)



Haier Biomedical  
International



Haier Biomedical  
International



@haiermedicalint



Haier Biomedical  
International



Haier Biomedical  
International

## Specifications



### Technical Data

Cabinet Type	Upright, Double Solid Door
Climate Class	N
Cooling Type	Forced Air Cooling
Defrost Mode	Auto
Refrigerant	HC
Noise (dB (A))	43

### Performance

Temp Range (°C)	2-8
-----------------	-----

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply (V/Hz)	220~240/50/60
Electrical Current (A)	2.7

### Dimensions

Capacity (L/Cu.Ft)	890/31.4
Net/Gross Weight (approx)	207/259(kg) 456.4/571.1(lbs)
Interior Dimensions (W*D*H)	1030*590*1425 (mm) 40.6*23.2*56.1 (in)
Exterior Dimensions (W*D*H)	1130*755*1980 (mm) 44.5*29.7*78.0 (in)
Packing Dimensions (W*D*H)	1250*850*2170 (mm) 49.2*33.5*85.4 (in)
Container Load (20'/40'/40'H)	10/22/22



### Alarms

High/Low Temp	Y
Remote Alarm	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

### Accessories

Caster	Y
Foot	Y
Porthole	Y
Shelves/Basket	12/-
Temperature Recorder	N/A
USB Interface	Optional
Printer	Optional

### Others

Certificate	CE
-------------	----