


# -40°C Biomedical Freezer



**DW-40L508**

## Scope of Application:

Designed to store vaccines, blood plasma and many other biological materials. Installations can be found in research institutions and clinical sites in the life science, pharma, biotech, medical and electronics sectors.

## Innovative Design

- ◉ Durable Refrigeration
- ◉ Ergonomic Design
- ◉ Environmentally-friendly
- ◉ Multiple Alarm Functions
- ◉ Large Capacity with Excellent Temperature Uniformity

### Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,  
Qingdao, 266109, P.R. China  
E-mail: [inquiry@haierbiomedical.com](mailto:inquiry@haierbiomedical.com)  
Website: [www.haiermedical.com](http://www.haiermedical.com)



Haier Biomedical  
International



Haier Biomedical  
International



@haiermedicalint



Haier Biomedical  
International



Haier Biomedical  
International

## Product Features



### Durable Refrigeration

Shelf style evaporator, fast cooling speed



### Six Alarm Functions

High/low temperature, power failure, high ambient temp, sensor error, low battery, remote alarm



### Drawer Structure

Multilayer drawer, easy to store and sort

## Product Advantages



### Lockable Design



### Hydrocarbon Environmentally Friendly



### High Efficiency and Energy Saving



### Test Hole Design



### Foot/Caster Design



### Intelligent Alarm

## Specifications

Model		DW-40L508		
Technical Data	Cabinet Type	Upright, Double Door		
	Climate Class	SN N		
	Cooling Type	Direct Cooling		
	Defrost Mode	Manual		
	Refrigerant	HC	CFC-Free	
	Sound Level (dB(A))	41	45	
Performance	Cooling performance (°C)	-40		
	Temperature Range (°C)	-20~-40		
Control	Controller	Microprocessor		
	Display	LED		
Electrical Data	Power Supply (V/Hz)	220-240/50	208-230/60	
	Electrical Current (A)	2.9	5.6	
Dimensions	Capacity (L/Cu.Ft)	490/17.3		
	Net/Gross Weight (approx)	kg	164/200	
		lbs	361.6/440.9	
	Interior Dimension (W*D*H)	mm	685*610*1228	
		in	27.0*24.0*48.3	
	Exterior Dimension (W*D*H)	mm	940*845*1860	
		in	37.0*33.3*73.2	
	Packing Dimension (W*D*H)	mm	1005*940*2040	
in		39.6*37.0*80.3		
Alarms	Container Load (20'/40'/40'H)	10/22/24		
	High/Low Temperature	Y		
	Remote Alarm	Y		
	Power Failure	Y		
	Sensor Error	Y		
	Low Battery	Y		
	High Ambient Temperature	Y		
	Dirty Condenser	N		
	Door Ajar	N		
	Accessories	Caster	Y	
		Leveling Feet	Y	
Porthole		Y		
Drawers/Inner Doors		10/N		
Adjustable Shelves		N		
USB Interface		Optional		
RS485		Optional		
NFC		N		
Others	Temperature Recorder	Optional		
	Certification	CE		

\*Product appearance and specifications are subject to change without notice.