

Biosafety Transport Container



FS-5 12 18 35 54LS

Scope of Application :

Applicable for cold chain transportation of samples, reagents, blood products, specimens and drugs, especially for the safe transportation of Category A (UN2814, UN2900), Category B (UN3373) infectious substances as well as PCR laboratory testing.

Innovative Design

- Safe, Secure and Reliable
- User-friendly, Easy to Operate
- Real-time Temperature Display
- Repeat Usage, Time and Cost Savings

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint





Haier Biomedical
International





Haier Biomedical
International


Product Features


- 


The sterilization process is easy, and the container can be used repeatedly
- 


High temperature outer case deformity protection for 10 years
- 

Manufactured with high-density polyethene and polypropylene material, which is non-toxic, UV resistant, non-fading and safe for direct food contact
- 

The container is filled with non-fluorinated polyurethane foam insulation material without any chemical materials, which is eco-friendly and without odor
- 

Easy to use: Smart-designed ice box, fixed clip, storage partition
- 

Integrated handles on both sides of the container, convenient, firm and durable
- 

Professional sealing with superior heat insulation and preservation effect
- 

Optional whole process temperature monitoring, which makes transportation more secure

Technical Parameters

Model	FS-5LS	FS-12LS	FS-18LS	FS-35LS	FS-54LS
QTY of Sealed Tanks (Pcs)	1	2	3	5	16
QTY of Sampling Tubes (Pcs)	20	40	60	100	320
Temperature Range (°C)	0-10	0-10	0-10	0-10	0-10
Duration (h)	≥6	≥24	≥24	≥24	≥24
Exterior Dimensions (W*D*H)(mm)	340*234*275	440*262*308	440*262*308	570*305*365	615*395*450
Interior Dimensions (W*D*H)(mm)	235*145*170	325*170*225	363*203*244	475*245*270	520*310*355
Net Weight without Tanks (kg)	4.8	6.2	6.5	9.5	18.5