



HYC-290

Pharmacy Refrigerator

Scope of Application

This product is suitable for pharmacy, pharmaceutical manufacturers, anti-epidemic stations, health centers and hospitals as key examples, which is used to store biological products and articles that need to be stored at 2~8°C.

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International

Product Features



Frequency Conversion Energy Saving

Variable frequency compressor with optimized cooling system, more energy efficient.



Environmentally Friendly

Hydrocarbon compressor and refrigerant, the refrigerant does not contain any fluorine or chlorine which is super friendly to the environment.



Large Space, Superior Uniformity

The space is efficiently used, and the utilization rate is high and the temperature uniformity is $\pm 3^{\circ}\text{C}$.



Low Noise

Optimized system and noise reduction design to reduce the noise level to 43dB(A) (220V/50Hz).



Stable performance, multiple alarm modes

Three alarm modes (light, sound, remote); six types of alarm.



Microprocessor Control System



Self Closing Door with 90° Stay Open Feature



Precise Temperature Control



Security lock



Optional USB and Printer

Specifications

Model		HYC-290	
Technical Data	Cabinet Type	Upright, Single Door	
	Climate Class	ST	
	Cooling Type	Forced Air Cooling	
	Defrost Mode	Auto	
	Refrigerant	HC	
	Sound Level (dB(A))	43	
Performance	Temperature Range ($^{\circ}\text{C}$)	2-8	
Control	Controller	Microprocessor	
	Display	LED	
Electrical Data	Power Supply (V/Hz)	115/60	220-240/50/60
	Electrical Current (A)	2.2	1.8
Dimensions	Capacity (L/Cu.Ft)	290/10.2	
	Net/Gross Weight (approx)	kg	105/127
		lbs	231.6/280.1
	Interior Dimensions (W*D*H)	mm	530*555*1080
		in	20.9*21.9*42.5
	Exterior Dimensions (W*D*H)	mm	665*710*1665
		in	26.2*28.0*65.6
	Packing Dimensions (W*D*H)	mm	700*810*1810
in		27.6*31.9*71.3	
Container Load (20'/40'/40'H)		22/44/44	
Alarms	High/Low Temperature	Y	
	Remote Alarm	Y	
	Power Failure	Y	
	Sensor Error	Y	
	Low Battery	Y	
	Door Ajar	Y	
Accessories	Caster	Y	
	Foot	Y	
	Porthole	Y	
	Shelves/Drawers	5/-	
	USB Interface	Y	Optional
	Printer	Optional	
Others	Certification	UL	CE