



HYC-940

Pharmacy Refrigerator



Six Alarm Functions

High/low temperature, remote alarm, power failure, sensor error, low battery, door ajar



Air Cooling, Auto Defrost

Better cabinet temperature uniformity, auto defrost



Adjustable Shelves

The shelves are adjustable to any height to fulfill different requirements



Heated Double Pane Glass Door

Clear display of stored products without condensation at 32°C and 85% relative humidity condition



Self Closing Door

Energy conserving, maintains temperature uniformity to protect valuable samples

Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International

Specifications

Model		HYC-940		
Technical Data	Cabinet Type	Upright, Double Glass Door		
	Climate Class	N		
	Cooling Type	Forced Air Cooling		
	Defrost Mode	Auto		
	Refrigerant	R290		
	Sound Level (dB(A))	43	45	
Performance	Temperature Range (°C)	2~8		
Control	Controller	Microprocessor		
	Display	LED		
Electrical Data	Power Supply (V/Hz)	220~240/50/60	115/60	
	Electrical Current (A)	3.5	7.5	
Dimensions	Capacity (L/Cu.Ft)		890/31.4	
	Net/Gross Weight (approx)	kg	227/279	
		lbs	500.4/615.1	
	Interior Dimensions (W*D*H)	mm	1030*590*1425	
		in	40.6*23.2*56.1	
	Exterior Dimensions (W*D*H)	mm	1130*755*1980	
		in	44.5*29.7*78.0	
	Packing Dimensions (W*D*H)	mm	1250*850*2170	
in		49.2*33.5*85.4		
Container Load (20'/40'/40'H)		10/22/22		
Alarms	High/Low Temperature	Y		
	Remote Alarm	Y		
	Power Failure	Y		
	Sensor Error	Y		
	Low Battery	Y		
	Door Ajar	Y		
Accessories	Caster	Y		
	Foot	Y		
	Porthole	Y		
	Shelves/Drawers	12/-		
	USB Interface	Optional		
	Printer	Optional		
	Chart Recorder	N/A		
Others	Certification	CE	UL	

*Haier Biomedical reserves the right to change products and specifications without prior notice.