

# **Benchtop High-Speed Refrigerated Microcentrifuge** LX-165T2R



# **Scope of Application**

Mainly used for nucleic acid extraction, DNA/RNA separation, purification, virus separation and protein precipitation in the laboratory. Widely used in universities, third party testing organizations, disease prevention, medical and health, biological pharmaceutical and other related fields.



Short pre-cooling time of only 5 minutes from room temperature to 4°C



Precise temperature and centrifugal speed control ensures optimal separation results



Optimized sound design for low noise levels (≤46dB(A))

Qingdao Haier Biomedical Technology Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone, Qingdao, 266109, P.R. China E-mail: inquiry@haierbiomedical.com Website: www.haiermedical.com

















# Benchtop High-Speed Refrigerated Microcentrifuge 🔻





### **Product Features**



Short pre-cooling time of only 5minutes from room temperature to 4°C



Fast sample processing at maximum speeds of up to 16,500 rpm (25,872 x a)



Precise temperature and centrifugal speed control ensures optimal separation results



Use of brushless induction motors to provide maintenance-free operation, eliminating carbon particle emissions



Optimized sound design for low noise levels (≤46dB(A))



Automatic radio-frequency rotor identification (RFID) and rotor life monitoring for maximum operational safety



Real time loading imbalance detection, immediate shutdown protection in case of imbalance



### Intuitive Interface



The 7-inch 45° glove-friendly touchscreen, the wide viewing angle more comfortable and practical for users of all heights



Further enhanced with access controls such as user login with password protection, the independent program data of each user is managed separately



Quick recall of up to 1000+ programs with alpha-numeric naming to minimize time between runs and for quick run start by routine users



The speed, centrifugal force, time and temperature can be changed by user during the centrifugation process



2 acceleration + 3 deceleration stages





## **Product Parts**



Intelligent lid locking monitoring and automatic self-locking to ensure safe operation



Viewing port in the lid for observing the operating condition of the rotor



Drainage tray to prevent internal water accumulation and improve corrosion resistance



to ensure safe use



316 stainless steel centrifuge chamber is rust-free and e asy to clean



Emergency lid release, pry out the white plastic plug on the side plate with a tool and pull the attached release cord at least 7 times to open the lid



USB port for exporting operation record including speed, time and temperature, operation curve with alarm messages etc.



# **Specifications**

Model	Maximum Capacity (ml)	Max Speed (rpm)	Max RCF(xg)	Speed Accuracy	Time Setting Range	Precooling Function	Temperature Setting Range(°C)	Display
LX-165T2R	12*10ml	16500	25872	±10rpm	1s~99h59m59s (1-second increment)	Yes	-20∼40, adjustable to 1	7 inch touch screen

Cooling	Noise Level (dB(A))	Power Supply (V/Hz)	Power(W)	Exterior Dimensions (W*D*H)(mm)	Net Weight (kg)	Packing Dimensions (W*D*H)(mm)	Gross Weight(kg)	Certifications
Direct Cooling	≤46	AC220/50/60	480	340*625*248	34	425*795*430	48	CE

# Benchtop High-Speed Refrigerated Microcentrifuge 🔻















BE No.	BE12C0Z00	BE11D0Z00	BE12D0Z00	BE12E0Z00	BE12DR000
Туре	Plastic cover	Sealed	Sealed	Sealed	Sealed
Max Capacity (ml)	24*1.5/2ml	24*1.5/2ml	12*5ml	4*8*0.2ml	12*10ml
Max RCF (xg)	25872	21382	21407	18891	2753
Max Speed (rpm)	16500	15000	15000	15000	5000
Tubes per Rotor	24	24	12	4	12
Tube				9999999	
Max Tube Capacity (ml)	1.5/2ml	1.5/2ml	5ml	0.2 ml 8-strip PCR	10ml
Max Tube Dimensions (Ø x L mm)	11*38	11*38	13.5*75	6*18	15*105
Tube Type	Round /Conical bottom	Round /Conical bottom	Round bottom	PCR	Large blood collection tube

### Adapter(Sets of 50) for BE12DR000

BE No.	BE12JC002	
Max Capacity (ml)	5ml	
Tube		
Max Tube Dimensions (Ø x L mm)	13*75	
Tube Type	Small blood collection tube	