

# Biosafety Transport Container



FS-12 18 35 54LS

## Scope of Application :

Applicable for cold chain transportation of samples, reagents, blood products, specimens and drugs, especially for the safe transportation of Category A (UN2814, UN2900), Category B (UN3373) infectious substances as well as PCR laboratory testing.

## Innovative Design

- Safe, Secure and Reliable
- Real-time Temperature Display
- User-friendly, Easy to Operate
- Repeat Usage, Time and Cost Savings

### Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,  
Qingdao, 266109, P.R. China  
E-mail: [inquiry@haierbiomedical.com](mailto:inquiry@haierbiomedical.com)  
Website: [www.haiermedical.com](http://www.haiermedical.com)



Haier Biomedical  
International



Haier Biomedical  
International



@haiermedicalint





Haier Biomedical  
International





Haier Biomedical  
International


## Product Features


- 


The sterilization process is easy, and the container can be used repeatedly
- 


High temperature outer case deformity protection for 10 years
- 

Manufactured with high-density polyethene and polypropylene material, which is non-toxic, UV resistant, non-fading and safe for direct food contact
- 

The container is filled with non-fluorinated polyurethane foam insulation material without any chemical materials, which is eco-friendly and without odor
- 

Easy to use: Smart-designed ice box, fixed clip, storage partition
- 

Integrated handles on both sides of the container, convenient, firm and durable
- 

Professional sealing with superior heat insulation and preservation effect
- 

Optional whole process temperature monitoring, which makes transportation more secure

## Technical Parameters

| Model                             | FS-12LS     | FS-18LS     | FS-35LS     | FS-54LS     |
|-----------------------------------|-------------|-------------|-------------|-------------|
| Temperature Range(°C)             | 2-10°C      | 2-10°C      | 2-10°C      | 2-10°C      |
| Duration(H)                       | ≥24h        | ≥24h        | ≥24h        | ≥24h        |
| Exterior Dimensions (W*D*H)(mm)   | 425*262*305 | 425*262*305 | 570*305*365 | 615*415*450 |
| Interior Dimensions (W*D*H)(mm)   | 325*170*225 | 375*215*250 | 475*245*270 | 510*315*355 |
| QTY of Sealed Tanks(Pcs)          | 2           | 3           | 5           | 14          |
| QTY of Samipling Tubes(Pcs)       | 40          | 60          | 100         | 280         |
| Qty of Icepacks(Pcs)              | 4           | 4           | 8           | 3           |
| Net Weight (include icepacks)(Kg) | 5.2         | 6.5         | 9           | 20          |
| Certification                     | CE          | CE          | CE          | CE          |