

Medical Transport Cooler



FS-12LS/18LS/35LS/54LS/
100L/110L

Scope of Application :

Cold chain transportation for samples, reagents, blood products, specimens and drugs.

Innovative Design

- ◉ Lightweight design, easy to handle
- ◉ Real-time display of inside temperature, more intuitive
- ◉ Non-pollution material with superior temperature control performance

Qingdao Haier Biomedical Co., Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266111, P.R. China
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com
Dec, 2025



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International

Product Advantages



High strength material, anti-collision, non-toxic and harmless, UV resistant.



The outer-unit body is equipped with an integrated handle, which is firmly installed and durable.



Real-time display of temperature inside the cooler, making the temperature monitoring more intuitive.



Specially positioned ice packs, partitioned storage, easy to operate.



Professional & innovative sealed design ensuring zero cold leakage to achieve safe, secure and reliable long-term cooling preservation.

Parameters

Model	FS-12LS	FS-18LS	FS-35LS	FS-54LS	FS-100L	FS-110L
Temperature Range (°C)	0-10	0-10	0-10	0-10	0-10	0-10
Duration (h)	≥24	≥24	≥24	≥24	≥24	≥24
Exterior Dimensions (W*D*H)(mm)	425*262*305	425*262*305	570*305*365	615*415*450	800*440*440	870*520*400
Interior Dimensions (W*D*H)(mm)	325*170*225	363*203*244	475*245*270	510*315*355	700*385*340	790*460*330
Qty of Ice packs (Pcs)	4	4	4	3	6	8
Net Weight (include ice packs)(kg)	4.8	5.2	7.5	15	19.6	24
Certification	CE	CE	CE	CE	CE	CE

Note: -25~-15°C and 15~25°C are available with optional ice packs. The standard ice packs can only achieve 0~10°C temperature storage.