

Disposable Plasma Apheresis Set

Scope of Application:

This product is used for an intermittent flow centrifugal plasma separator for the separation, collection, or removal of human plasma.

Product Composition:

centrifuge bowl, plasma tube, plasma bag, venous fistula needle set, sterilized with ethylene oxide for single use.

Advantages:

The plasma separation is fast, efficient and low noise, which can be used in the application scenario with high specification requirements. It is composed of several relatively independent parts, modular design, easy operation and high reliability; Each component module can be freely combined and used in combination to meet different usage scenarios and individual needs.



Approval Number: 20150022

Qingdao Haier Biomedical Co., Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266111, P.R. China
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com
Jan, 2026



Haier Biomedical



Haier Biomedical
International



Haier Biomedical
International



Haier Biomedical
International




Haier Biomedical
International



Centrifuge Bowl

Low noise;
High light transmittance material;
Few red blood cells adhere to the wall;
Patented technology for preventing temperature rise.



Tube

High-quality filters;
PVC main material;
Multiple specifications are available;
High light transmittance and high flexibility.



Plasma Bag

Mature plastic sealing process;
Multiple specifications are available;
It can be used alone as a liquid storage bag;
High light transmittance and high flexibility PVC main material.



Venous Fistula Needle Set

Multiple specifications are available;
It can be used as a blood collection needle alone;
High-quality imported, grind-and-molded stainless-steel needle;
High light transmittance and high flexibility PVC main material.