



Glass and Solid Door
Models Available

HYC-80GE



Desk/under-desk Placement,
Reversible Door



Auto-defrost, Frost-free



One-touch High/
Low Temp Check



Manual Reset



Easy Operation

Pharmacy Refrigerator

Scope of Application

Used for the storage of medicines, vaccines, reagents, biological products, and other items that require a temperature between 2–8°C in university laboratories, healthcare institutions, pharmaceutical companies, nursing homes, medical aesthetics clinics, etc.

Innovative Design Features

- Left/Right Reversible Door
- Low-noise Operation
- Low Energy Consumption

Qingdao Haier Biomedical Co., Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266111, P.R. China
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com
Feb, 2026



Haier Biomedical



Haier Biomedical
International



Haier Biomedical
International



Haier Biomedical
International



Haier Biomedical
International

Product Features



85L Net Volume, Reversible Door Design



All-round Safety Alarm System, Built-in Remote Alarm Interface



36dB(A) Ultra-quiet, Air-cooled Design, Fluorine-free Inverter Compressor



One-click Query for High & Low Temperatures, Manual Reset; Simple and Easy to Operate



CE Certified; DIN DIN Compliant




Humanized Design: Anti-misoperation Touch Power Switch, Rear-mounted, Optional USB, Auto Defrost

Product Details



Product Dimension



Specifications

	Model		HYC-80GE
Technical Data	Cabinet Type		Counter-top/Undercounter Glass Door
	Climate Class		N, 16-32°C
	Cooling Type		Forced Air Cooling
	Defrost Mode		Auto
	Refrigerant		R600a
	GWP (Global Warming Potential)		20
	Sound Level (dB(A))		42
Performance	Temperature Range (°C)		2-8
Control	Will the fan stop automatically when opening the door		Y
	Controller		Microprocessor
	Display		LED
Electrical Data	Power Supply (V/Hz)		220-240/50/60
	Power (W)		180
Dimensions	Capacity (L/Cu.Ft)		85 / 3.0
	Net/Gross Weight (approx)	kg	/
		lbs	/
	Interior Dimensions (W*D*H)	mm	405*405*555
		in	15.9*15.9*21.9
	Exterior Dimensions (W*D*H)	mm	495*550*740
		in	19.5*21.7*29.1
	Packing Dimensions (W*D*H)	mm	590*650*830
		in	23.2*25.6*32.7
Alarms	Container Load (20'/40'/40'H)		54/120/180
	High/Low Temperature		Y
	High Ambient Temperature		Y
	Door Ajar		Y
	Power Failure		Y
	Low Battery		Y
	Sensor Error		Y
	USB Communication Failure		Optional
	Remote Alarm		Y
	Audible and Visual Alarm		Y
Accessories	Roller		2
	Foot		2
	Porthole		1 R
	Shelves/Basket		2/1
	Printer		N/A
	Chart Recorder		N/A
	USB Interface		Optional
	RS485		Optional
Others	Certification		CE, DIN13277 (Optional)

*Haier Biomedical reserves the right to change products and specifications without prior notice.